

Human IL-6RA/CD126 protein (Recombinant) (C-hFc & His) (STJP002757)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	3570
Gene Symbol	IL6R
UniProt ID	IL6RA_HUMAN
Immunogen Region	Leu20-Asp358
Immunogen Sequence	LAPRRCPAQEVARGVLTSLP GDSVTLTCPGVEPEDNATVH WVLKPAAGSHPSRWAGMGR RLLRSVQLHDSGNYSYRA GRPAGTVHLLVDVPPEEPQL SCFRKSPLSNVCEWGPRST PSLTTKAVLLVRKFQNSPAE DFQEPQYSQESQKFSCQLA VPEGDSSFYIVSMCVASSVG SKFSKTQTFQGCGILQPDPP ANITVTAVARNPRWLSVTWQ DPHSWNSSFYRLRFELRYR

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Recombinant Human IL-6RA/CD126 Protein was determined by SDS-PAGE with Coomassie Blue, showing a band at 80-100kDa.



Immobilized recombinant Human IL6 at 1 Mu g/mL (100 Mu L/well) can bind Human IL6RA with a linear range of 1-5.4 ng/mL.