

Human Ephrin-A1/EFNA1 protein (Recombinant) (C-His) (STJP002759)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES


Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	1942
Gene Symbol	EFNA1
UniProt ID	EFNA1_HUMAN
Immunogen Region	Asp19-Ser182
Immunogen Sequence	DRHTVFWNSSNPKFRNEDYT IHVQLNDYVDIICPHYEDHS VADAAMEQYILYLVEHEEYQ LCQPQSKDQVRWQCNRPSAK HGPEKLSEKFQRFPTLGK EFKEGHSYYYISKPIHQHED RCLRLKVTVSGKITHSPQAH DNPQEKRLAADDPEVRLHS IGHS

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



Active Recombinant Human Ephrin-A1/EFNA1 Protein was determined by SDS-PAGE with Coomassie Blue, showing a band at 26kDa.



Immobilized Human EFNA1 at 1 Mu g/mL (100 Mu L/well) can bind Human Ephrin A1 Rabbit mAb with a linear range of 0. 098-2. 066 ng/mL.