

Human GITR Ligand/TNFSF18 protein (Recombinant) (C-hFc) (STJP002766)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	8995
Gene Symbol	TNFSF18
UniProt ID	TNF18_HUMAN
Immunogen Region	Gln50-Ser177
Immunogen Sequence	QLETAKEPCMAKFGPLPSKW QMASSEPPCVNKVSDWKLEI LQNGLYLIYGQVAPNANYND VAPFEVRLYKNKDMIQTLTN KSKIQNVGGTYELHVGDTID LIFNSEHQVLKNNTYWGII LANPQFIS

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---