

Mouse FGFR-4/CD334 protein (Recombinant) (C-hFc) (STJP002767)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES


Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	14186
Gene Symbol	Fgfr4
UniProt ID	FGFR4_MOUSE
Immunogen Region	Leu17-Asp366
Immunogen Sequence	LSLEASEEMEQEPCLAPILE QQEQVLTVALGQPVRLLCCGR TERGRHWYKEGSRLASAGRV RGWRGRLEIASFLPEDAGRY LCLARGSMTVVHNLTLMDL SLTSISNDEDPKTLSSSSSG HVYPQQAPYWTHPQRMEKKL HAVPAGNTVKFRCPAAGNPM PTIHWLKDQGAFFHGENRIGG IRLRHQHWSLVMESVVPDR GTYTCLVENSLSIRYSYLL DVLESPHRPILQAGLPAN

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Recombinant Mouse FGFR-4/CD334 Protein was determined by SDS-PAGE with Coomassie Blue, showing a band at 80-95kDa.



Immobilized Human FGF2 Protein at 5 Mu g/mL (100 Mu L/well) can bind Mouse FGFR-4/CD334 with a linear range of 0. 12-53. 83 ng/mL.