

Mouse CD200R1 protein (Recombinant) (C-hFc & Avi) (STJP002768)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES


Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	57781
Gene Symbol	Cd200r1
UniProt ID	MO2R1_MOUSE
Immunogen Region	Thr26-Pro238
Immunogen Sequence	TDKNQTTQNNSSSPLTQVNT TVSVQIGTKALLCCFSIPT KAVLITWIKLRGLPSCITIA YKVDTKTNETSCLGRNITWA STPDHSPQLQISAVTLQHEG TYTCETVTPEGNFEKNYDLQ VLVPPEVTYFPEKNRSVCE AMAGKPAAQISWSPDGDVCT TSESHSNGT VTVRSTCHWEQ NNVSDVSCIVSHLTGNQSL SIELSRGGNQSLRP

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Recombinant Mouse CD200R1 Protein was determined by SDS-PAGE with Coomassie Blue, showing a band at 75-90kDa.



Immobilized Human CD200 Protein at 2 Mu g/mL (100 Mu L/well) can bind mouse CD200 R1 with a linear range of 3. 9-92. 9ng/mL.