

## Human gp140/CDCP1/CD318 protein (Recombinant) (C-hFc) (STJP002770)

### GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

### PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">64866</a>
Gene Symbol	<a href="#">CDCP1</a>
UniProt ID	<a href="#">CDCP1_HUMAN</a>
Immunogen Region	Phe30-Leu666
Immunogen Sequence	FEIALPRESNITVLIKLGTP TLLAKPCYIVISKRHITMLS IKSGERIVFTFSCQSPENHF VIEIQKNIDCMGSPCFGEV QLQPSTSLPLNRTFIWDV KAHKSIGLELQFSIPRLRQI GPGESCPDGVTHSISGRIDA TVVRIGTFCSNGTVSRIKMQ EGVKMALHLPWFHPRNVSGF SIANRSSIKRLCIIESVFEG EGSATLMSANYPEGFPEDEL MTWQFVVP AHLRASVSFLN

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	-------------------------------------------------------------------------------------------------------------------------