

Mouse HMGB1 protein (Recombinant) (C-His) (STJP002772)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 M µm filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	15289
Gene Symbol	Hmgb1
UniProt ID	HMGB1_MOUSE
Immunogen Region	Met1-Glu215
Immunogen Sequence	MGKGDPPKPRGKMSSYAFFV QTCREEHKKKHPDASVNFSE FSKKCSERWKTMSAKEKGKF EDMAKADKARYEREMKTYIP PKGETKKKFKDPNAPKRPPS AFFLFCSEYRPKIKGEHPGL SIGDVAKKLGEMWNNTAADD KQPYEKKAALKKEYEKDIA AYRAKGKPDAAKKGVVKA EK SKKKKEEEDDEEEDDEEEEEE EEEEDDEEEDDDDE

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---