

Human IL-1 Rrp2/IL-1 R6 protein (Recombinant) (C-hFc) (STJP002773)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	8808
Gene Symbol	IL1RL2
UniProt ID	ILRL2_HUMAN
Immunogen Region	Asp20-Arg335
Immunogen Sequence	DGCKDIFMKNEILSASQPFA FNCTFPPITSGEVSVTWYKN SSKIPVSKIIQSRIHQDETW ILFLPMEWGD SGVYQCVIKG RDSCHRIHVNLTVFEKHWCD TSIGGLPNLSDEYKQILHLG KDDSLTCHLHFPKSCVLGPI KWKYKDCNEIKGERFTVLETR LLVSNVSAEDRGNYACQAIL THSGKQYEVNLNGITVSITER AGYGGSVPKIIPKNHSIEV QLGTTLIVDCNVTDTKDNT

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---