

Mouse CD83 protein (Recombinant) (C-His) (STJP002780)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 M µm filtered solution of 20 mM tris, 200 mM NaCl, pH 8.0
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	12522
Gene Symbol	Cd83
UniProt ID	CD83_MOUSE
Immunogen Region	Met22-Ala134
Immunogen Sequence	MREVTVACSETADLPCTAPW DPQLSYAVSWAKVSESGTES VEPESKQNSSFAPRRRAY SLTIQNTTICSSGTYRCALQ ELGGQRNLSGTVVLKVTGCP KEATESTFRKYRA

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---