

Mouse TREM-2 protein (Recombinant) (C-His) (STJP002784)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES


Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	83433
Gene Symbol	Trem2
UniProt ID	TREM2_MOUSE
Immunogen Region	Leu19-Ser171
Immunogen Sequence	LNTTVLQGMAGQSLRVSCY DALKHWGRRKAWCRQLGEEG PCQRRVSTHGVWLLAFLKKR NGSTVIADDTLAGTVITLK NLQAGDAGLYQCQSLRGREA EVLQKVLVEVLEDPLDDQDA GDLWVPEESSSEGAQVEHS TSRNQETSFPTTS

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



Recombinant Mouse TREM-2 Protein was determined by SDS-PAGE with Coomassie Blue, showing a band at 25-40 kDa.



Immobilized Mouse TREM-2 Protein at 1 Mu g/mL (100 Mu L/well) can bind TREM-2 Rabbit pAb with a linear range of 0.06-11.7 ng/mL.