

Mouse TNFRSF9/4-1BB/CD137 protein (Recombinant) (C-hFc) (STJP002786)

GENERAL INFORMATION	
Product Type	Proteins
Host / Source	HEK293 cells
PRODUCT PROPERTIES	
Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
TARGET INFORMATION	
Gene ID	21942
Gene Symbol	Tnfrsf9
UniProt ID	TNR9_MOUSE
Immunogen Region	Val24-Leu211
Immunogen Sequence	VQNSCDNCQPGTFCRKYNPV CKSCPPSTFSSIGGQPNCNI CRVCAGYFRFKFCSSTHNA ECECIEGFHCLGPQCTRCEK DCRPGQELTKQGCKTCSLGT FNDQNGTGVCRPWTNCSLDG RSVLKTGTTEKDVVCGPPVV SFSPSTTISVTPEGGPAFKK TTGAAQEEDACSCRCPQEEE GGGGGYEL
ADDITIONAL INFORMATION	
Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.