

## Mouse Sonic hedgehog N-product/SHH protein (Recombinant) (C-His) (STJP002787)

### GENERAL INFORMATION

|               |              |
|---------------|--------------|
| Product Type  | Proteins     |
| Host / Source | HEK293 cells |

### PRODUCT PROPERTIES

|                     |  |
|---------------------|--|
| Formulation         | Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.                                 |
| Storage Instruction | Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles. |

### TARGET INFORMATION

|                    |  |
|--------------------|--|
| Gene ID            | <a href="#">20423</a>  |
| Gene Symbol        | <a href="#">Shh</a>  |
| UniProt ID         | <a href="#">SHH_MOUSE</a>  |
| Immunogen Region   | Ala24-Gly198   |
| Immunogen Sequence | ACGPGRGFGKRRHPKLTPL AYKQFIPNVAEKT LGASGRY EGKITRNSERFKELTPNYNP<br>DIIFKDEENTGADRLMTQRC KDKLNALAISVMNQWPGVKL RVTEGWDEDGHHSEESLHYE<br>GRAVDITTSRDRRSKYGMLA RLAVEAGFDWVYYESKAHIIH CSVKAENSVAAKSGG |

### ADDITIONAL INFORMATION

|      |   |
|------|---|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|---|