

Mouse ALCAM/CD166 protein (Recombinant) (C-His) (STJP002796)

GENERAL INFORMATION	
Product Type	Proteins
Host / Source	HEK293 cells
PRODUCT PROPERTIES	
Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
TARGET INFORMATION	
Gene ID	11658
Gene Symbol	Alcam
UniProt ID	CD166_MOUSE
Immunogen Region	Trp28-Lys527
Immunogen Sequence	WYTVNSAYGDTIVMPCRLDV PQNLMFGKWKEYEKPDGSPVF IAFRSSTKKSVQYDDVPEYK DRLSLSENYTLSIANAKISD EKRFVCMVLTEDNVFEAPTL VKVFKQPSKPEIVNKAPFLE TDQLKKLGDCISRDSYPDGN ITWYRNGKVLQPVEGEVAIL FKKEIDPGTQLYTVTSSEY KTTRSDIQMPFTCSVTYYGP SGQKTIYSEQEIFDIYYPTQVTIQVLPKNAIKEGDN
ADDITIONAL INFORMATION	
Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.