

## Mouse Asparaginyl endopeptidase/Legumain protein (Recombinant) (C-His) (STJP002798)

### GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

### PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">19141</a>
Gene Symbol	<a href="#">Lgmn</a>
UniProt ID	<a href="#">LGMN_MOUSE</a>
Immunogen Region	Val18-Tyr435
Immunogen Sequence	VPVGDDPEDGGKHWVIVA GSNGWYNYRHQADACHAYQI IHRNGIPDEQIIVMMYDDIA NSEENPTPGVINRPNGTDV YKGVLDYTGEDVTPENFLA VLRGDAEAVKGGKSGKVLKS GPRDHVFIYFDHGATGILV FPNDDLHVKDLNKTIRYMYE HKMYQKMVFYIEACESGSM NHLPDDINVYATTAANPKES SYACYDEERGTYLGDWYSV NWMEDSDVEDLTKETLHKQ

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---