

Mouse SIRP-alpha/CD172a protein (Recombinant) (C-His) (STJP002799)

GENERAL INFORMATION	
Product Type	Proteins
Host / Source	HEK293 cells
PRODUCT PROPERTIES	
Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
TARGET INFORMATION	
Gene ID	19261
Gene Symbol	Sirpa
UniProt ID	SHPS1_MOUSE
Immunogen Region	Lys32-Asn373
Immunogen Sequence	KELKVTQPEKSVSVAAGDST VLNCTLTSLLPVGPIRWYRG VGPSRLLIYSFAGEYVPRIR NVSDTTKRNNMDFSIRISNV TPADAGIYYCVKFQKGSSEP DTEIQSGGGTEVYVLAKPSP PEVSGPADRGIPDQKVNFTC KSHGFSPRNITLKWFKDGQE LHPLETTVNP SGKNVSYNIS STVRVVLNSMDVNSKVICEV AHITLDRSPLRGIANLSNFI RVSP TVKVTQQSPTSMNQV
ADDITIONAL INFORMATION	
Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.