

Mouse Nectin-2/PVRL2/CD112 protein (Recombinant) (C-His) (STJP002800)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	19294
Gene Symbol	Nectin2
UniProt ID	NECT2_MOUSE
Immunogen Region	Gln32-Gly351
Immunogen Sequence	QDVRVRLPEVRGRLGGTVE LPCHLLPPTTERVSQVTWQR LDGTVVAAFHPSFGVDFPNS QFSKDRLSFVRARPETNADL RDATAFRGLRVEDEGNYTC EFATFPNGTRRGVTWLRVIA QPENHAAEAQEVITIGPQSVAV ARCVSTGGRPPARITWISSL GGEAKDTQEPGIQAGTVTII SRYSLVPVGRADGVKTCRV EHESFEEPILLPVTLVSRYP PEVSISGYDDNWYLRSEA

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---