

Human Lysyl hydroxylase 2/PLOD2 protein (Recombinant) (C-His) (STJP002805)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Supplied as a 0.22 µg/ml filtered solution in PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5352
Gene Symbol	PLOD2
UniProt ID	PLOD2_HUMAN
Immunogen Region	ALa26-Pro737
Immunogen Sequence	ADSEKPSIPTDKLLVITVA TKESDGFHRFMQSAKYFNVT VKVLGQGEEWRGGDGINSIG GGQKVRMLKEVMEHYADQDD LVVMFTECFDVIFAGGPPEV LKKFQKANHKVVFAADGILW PDKRLADKYVHVHIGKRYLN SGGFIGYAPYVNRIVQQWNL QDNDDQLFYTKVYIDPLKR EAINITLDHKCKIFQTLNGA VDEVVLKFEN GKARAKNTFY ETLPVAINGNGPTKILLNY

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---