

Mouse Clusterin/Apo-J/CLU protein (Recombinant) (C-His) (STJP002815)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	12759
Gene Symbol	Clu
UniProt ID	CLUS_MOUSE
Immunogen Region	Gly21-Glu448
Immunogen Sequence	GEQEVSDNELQELSTQGSRY INKEIQNAVQGVKHIKTLIE KTNAERKSLLSLEEAKKKK EDALETRDSEMKLKAFPEV CNETMMALWEECKPCLKHTC MKFYARVCRSGSGLVGQQL EFLNQSSPFYFWMNGDRIDS LLESDRQQSQVLDAMQDSFA RASGIIDTLFQDRFFARELH DPHYFSPIGFPHKRPFLYP KSRLVRSLMSPSHYGPPSFH NMFQPPFEMIHQAAQAMDV

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---