

Human CCL5/RANTES protein (Recombinant) (C-His) (STJP002819)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	Pichia

PRODUCT PROPERTIES

Formulation	Supplied as a 0.22 µg/ml filtered solution in PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	6352
Gene Symbol	CCL5
UniProt ID	CCL5_HUMAN
Immunogen Region	Ser24-Ser91
Immunogen Sequence	SPYSSDTPCCFAYIARPLP RAHIKEYFYTSGKCSNPAVV FVTRKNRQVCANPEKKWVRE YINSLEMS

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---