

## Human CCL1/l-309 protein (Recombinant) (No tag) (STJP002821)

### GENERAL INFORMATION

Product Type	Proteins
Host / Source	Pichia

### PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">6346</a>
Gene Symbol	<a href="#">CCL1</a>
UniProt ID	<a href="#">CCL1_HUMAN</a>
Immunogen Region	Lys24-Lys96
Immunogen Sequence	KSMQVPFSRCCFSFAEQEIP LRAILCYRNTSSICSNEGLI FKLKRGKEACALDTVGWVQR HRKMLRHCPKSRK

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---