

Human CCL13/MCP-4 protein (Recombinant) (No tag) (STJP002822)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	Pichia

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	6357
Gene Symbol	CCL13
UniProt ID	CCL13_HUMAN
Immunogen Region	Gln24-Thr98
Immunogen Sequence	QPDALNVPSTCCFTFSSKKI SLQRLKSYVITTSRCPQKAV IFRTKLGKEICADPKEKVVQ NYMKHLGRKAHTLKT

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---