

Human Afamin/AFM protein (Recombinant) (C-His) (STJP002838)

GENERAL INFORMATION	
Product Type	Proteins
Host / Source	HEK293 cells
PRODUCT PROPERTIES	
Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
TARGET INFORMATION	
Gene ID	173
Gene Symbol	AFM
UniProt ID	AFAM_HUMAN
Immunogen Region	Met1-Asn599
Immunogen Sequence	MKLLKLTGFIFFLFLTESL TLPTQPRDIENFNSTQKFIE DNIEYITIIAFAQYVQEATF EEMEKLVKDMVEYKDRCMAD KTLPECSKLPNNVLQEKICA MEGLPQKHNFHSCCSKVDAQ RRLCFFYNKKSDVGFLPPFP TLDPEEKQAYESNRESLLN HFLYEVARRRNPFVFAPTLT VAVHFEEVAKSCCEEQNKVN CLQTRAIPTVQYLKAFSSYQ KHVCGALLKFGTKVVHFIY
ADDITIONAL INFORMATION	
Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.