

Human CXCL12/SDF-1 protein (Recombinant) (No-Tag) (STJP002843)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	Pichia

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	6387
Gene Symbol	CXCL12
UniProt ID	SDF1_HUMAN
Immunogen Region	Lys22-Lys89
Immunogen Sequence	KPVSLSYRCPCRFFESHVAR ANVKHLKILNTPNCALQIVA RLKNNNRQVCIDPKLKWIQE YLEKALNK

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---