

Human IL-10 R alpha/IL10RA/CD210 protein (Recombinant) (C-RaFC) (STJP002852)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES


Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	3587
Gene Symbol	IL10RA
UniProt ID	I10R1_HUMAN
Immunogen Region	His22-Asn235
Immunogen Sequence	HGTLPSPPSVWFEEFFHH ILHWTPIPNQSESTCYEVAL LRYGIESWNSISNCSQTLSTY DLTAVTLDLYHSNGYRARVR AVDGSRHSNWTVTNTRFSVD EVTLTVGSVNLEIHNGFILG KIQLPRPKMAPANDTYESIF SHFREYEIAIRKVPGNFTFT HKKVKHENFSLTSGEVGEF CVQVKPSVASRSNKGWWSKE ECISLTRQYFTVTN

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Active Recombinant Human IL-10 R alpha Protein was determined by SDS-PAGE with Coomassie Blue, showing a band at 60-80 kDa.



Immobilized Human IL10 (Catalog: STJP000869) at 2 Mu g/mL (100 Mu L/well) can bind Human IL10RA (Catalog: STJP002852) with a linear range of 0. 1-2. 9 ng/mL.