

Human CCL4/MIP-1Beta protein (Recombinant) (C-hFc) (STJP002858)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	388372/6351
Gene Symbol	CCL4
UniProt ID	CCL4_HUMAN
Immunogen Region	Ala24-Asn92
Immunogen Sequence	APMGSDPPTACCFSYTARKL PRNFVVDYYETSSLCSQPAV VFQTKRSKQVCADPSESWWQ EYVYDLELN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	-------------------------------------------------------------------------------------------------------------------------