

Rat VEGFR-2/KDR/CD309 protein (Recombinant) (C-His) (STJP002862)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	25589
Gene Symbol	Kdr
UniProt ID	VGFR2_RAT
Immunogen Region	Ala20-Glu760
Immunogen Sequence	ASVGLPGDSLHPPKLSTQKD ILTILANTTLQITCRGQRDL DWLWPNTPRDSEERVLVTEC GDSIFCKTLTVPRVVGNDTG AYKCFYRDTDVSSIVYVYVQ DHRSPFIASVSDEHGIVYIT ENKNKTVVIPCRGSISNLNV SLCARYPEKRFVDPGNRISW DSEKGFIPSYMISYAGMVF CEAKINDETYQSIMYIVLVV GYRIYDVVLSPPEHELSAG EKLVLNCTARTELVGLDF

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---