

Mouse LILRB4/CD85k/ILT3 protein (Recombinant) (C-6His) (STJP002878)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	14728
Gene Symbol	Lilrb4a
UniProt ID	LRB4A_MOUSE
Immunogen Region	Gly24-Lys238
Immunogen Sequence	GHLPKPIIWAEPGSVIAAYT SVITWCQGSWEAQYYHLYKE KSVNPWDTQVPLETRNKAKF NIPSMTTSYAGIYKCYYESA AGFSEHSDAMELVMTGAYEN PSLSVYPSSNVTSGVSI SFS CSSSIVFGRFILIQEGKHGL SWTLDSQHQANQPSYATFVL DAVTPNHNGTFR CYGYFRNE PQVWSKPSNSLDLMISETKD QSSTPTEDGLETYQK

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---