

Mouse DLL4 protein (Recombinant) (C-6His) (STJP002880)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	54485
Gene Symbol	DII4
UniProt ID	DLL4_MOUSE
Immunogen Region	Gly27-Ala532
Immunogen Sequence	GSGIFQLRLQEFVNQRGMLA NGQSCEPGCRTFFRICLKHF QATFSEGPCTFGNVSTPVLG TNSFVVRDKNSGSGRNPLQL PFNFTWPGTFSLNIQAWHTP GDDLRPETSPGNSLISQIII QGSLAVGKIWRDTEQNDTLT RLSYSYRVICSDNYYGESCS RLCKKRDDHFGHYECQPDGS LSCLPGWTGKYCDQPICLSG CHEQNGYCSKPDECICRPGW QGRLCNECIPHNGCRHGTC

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---