

Mouse CD276 protein (Recombinant) (C-6His) (STJP002881)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES


Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	102657
Gene Symbol	Cd276
UniProt ID	CD276_MOUSE
Immunogen Region	Val 29-Phe 244
Immunogen Sequence	VEVQVSEDPVVALVDTDATL RCSFSPEPGFSLAQLNLIWQ LTDTKQLVHSFTEGRDQGSA YSNRTALFPDLLVQGNASLR LQRVRVTDEGSYTCFVSIQD FDSAASVSLQVAAPYSKPSMT LEPNKDLRPGNMVTITCSSY QGYPEAEVFWKDGQGVPLTG NVTTSQMANERGLFDVHSVL RVVLGANGTYSCLVRNPVLQ QDAHGSVTITGQPLTF

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Recombinant Mouse CD276 Protein was determined by SDS-PAGE with Coomassie Blue, showing a band at 35-45 kDa.



Immobilized Mouse CD276 (Catalog: STJP002881) at 2 Mu g/mL (100 Mu L/well) can bind CD276 Rabbit pAb (Catalog: STJ231053) with a linear range of 1-274 ng/mL.