

Rat PCSK9 protein (Recombinant) (C-6His) (STJP002884)

GENERAL INFORMATION	
Product Type	Proteins
Host / Source	HEK293 cells
PRODUCT PROPERTIES	
Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
TARGET INFORMATION	
Gene ID	298296
Gene Symbol	Pcsk9
UniProt ID	PCSK9_RAT
Immunogen Region	Gln 31-Gln 691
Immunogen Sequence	QDEDGDYEELMLALPSQEDS LVDEASHVATATFRRCSKEA WRLPGTYVVVLMEETQRLQV EQTAHRLQTWAARRGYVIKV LHVFDLFPGLVKMSSDLL GLALKLPHVEYIEEDSLVFA QSIPWNLERIIPAWQQTEED SSPDGSSQVEVYLLDTSIQS GHREIEGRVTITDFNSVPEE DGTRFHRQASKCDSHGTHLA GVVSGRDAGVAKGTSLSLRL VLNCQKGKTVSGTLIGLEF
ADDITIONAL INFORMATION	
Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.