

Human Betacellulin/BTC protein (Recombinant) (C-hFc) (STJP002897)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	685
Gene Symbol	BTC
UniProt ID	BTC_HUMAN
Immunogen Region	Asp32-Tyr111
Immunogen Sequence	DGNSTRSPETNGLLCGDPEE NCAATTTQSKRKGHFRCPK QYKHYCIKGRCRFVVAEQTP SCVCDEGYIGARCERVDLFY

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--