

Human CXCL9/MIG protein (Recombinant) (No tag) (STJP002906)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	E.coli

PRODUCT PROPERTIES


Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	4283
Gene Symbol	CXCL9
UniProt ID	CXCL9_HUMAN
Immunogen Region	Thr23-Thr125
Immunogen Sequence	MTTASTSQVRQNYHQDSEAA INRQINLELYASYVYLSMSY YFDRDDVALKNFAKYFLHQS HEEREHAEKLMKLNQRGGR IFLQDIKKPCDDWESGLNA MECALHLEKNVNSLLELHK LATDKNDPHLCDFIETHYLN EQVKAIKELGDHVTNLRKMG APESGLAEYLFDKHTLGSD NES

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



Recombinant protein Human CXCL9 was determined by SDS-PAGE under reducing conditions with Coomassie Blue, showing a band at 15-25 kDa.