

Human IL-37/IL-1F7 protein (Recombinant) (No-Tag) (STJP002909)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	E.coli

PRODUCT PROPERTIES


Formulation	Lyophilised from a 0.2 µm filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	27178
Gene Symbol	IL37
UniProt ID	IL37_HUMAN
Immunogen Region	Lys27-Asp192
Immunogen Sequence	NPAAEAEFNILLATDSYKVT HYKQYPPNTSKVSYFECRE KKTENSKLRKVKYEETVFYG LQYILNKYLKGGVVTKEKIQ EAKDVYKEHFQDDVFNEKGW NYILEKYDGHLPKAVPE GFVIPRGNVLFVTVENTDPEC YWLTNWIIETLVQSWYPITV ATNSREQKKILAKYLLETSG NLDGLEYLKLDHFGYRGVSSQ ETAGIGASAHLVNFKGTDTV AGLALIKKYYGKDPVPGY

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



Recombinant protein Human Interleukin-37/IL-37 was determined by SDS-PAGE under reducing conditions with Coomassie Blue, showing a band at 18 kDa.