

## Human FABP3/H-FABP protein (Recombinant) (N-His & No-Tag) (STJP002911)

### GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

### PRODUCT PROPERTIES


Formulation	Lyophilised from 0.22 Mu m filtered solution in PBS (pH 7.4). Normally 8% Trehalose is added as protectant before lyophilisation.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">2170</a>
Gene Symbol	<a href="#">FABP3</a>
UniProt ID	<a href="#">FABPH_HUMAN</a>
Immunogen Region	Val2-Ala133
Immunogen Sequence	PRPPGRPQLSLQELRREFTV SLHLARKLLSEVRGQHRFA ESHLPGVNLYLLPLGEQLPD VSLTFQAWRRLSDPERLCFI STTLQPFHALLGGLGTQGRW TNMERMQLWAMRLDLRDLQR HLRFQVLAAGFNLPPEEEEE EEEEEERKGLLPGALGSAL QGPAQVSWPQLLSTYRLLHS LELVLSRAVRELLLSKAGH SVWPLGFPTLSPQP&RKGPP AALTLPVQCQRASRYPIAV

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



Recombinant Human FABP3/H-FABP Protein was determined by SDS-PAGE with Coomassie Blue, showing a band at 15-20 kDa.