

Mouse Alpha-Synuclein/SNCA protein (Recombinant) (C-His) (STJP002912)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES


Formulation	Lyophilised from a 0.2 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	20617
Gene Symbol	Snca
UniProt ID	SYUA_MOUSE
Immunogen Region	Met1-Ala140
Immunogen Sequence	VWEKTVNTEENVYATLGSDV NLTCQTQTVGFFVQMOWSKV TNKIDLIAYVHPQYGFYCA YGRPCESLVTFETEPENGSKW TLHLRNMSCSVSGRYECMLV LYPEGIQTKIYNLLIQTHVT ADEWNSNHTIEIEINQTLIEI PCFQNSSSKISSEFTYAWSV EDNGTQETLISQNHLSNST LLKDRVKLGTDYRLHLSPVQ IFDDGRKFSCHIRVGNPKIL RSSTTVKVFAPKEIPVIVE

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



Recombinant Mouse SNCA/Alpha-Synuclein Protein was determined by SDS-PAGE with Coomassie Blue, showing a band at 15-25 kDa.