

## Human S100-B protein (Recombinant) (N-hFc) (STJP002915)

### GENERAL INFORMATION

|               |              |
|---------------|--------------|
| Product Type  | Proteins     |
| Host / Source | HEK293 cells |

### PRODUCT PROPERTIES


|                     |                                                                                                |
|---------------------|------------------------------------------------------------------------------------------------|
| Formulation         | Lyophilised from a 0.2 µm filtered solution of 20mM PBS, 150mM NaCl, pH 7.4.                   |
| Storage Instruction | Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles. |

### TARGET INFORMATION

|                    |                                                                                                                                                                         |
|--------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Gene ID            | <a href="#">6285</a>                                                                                                                                                    |
| Gene Symbol        | <a href="#">S100B</a>                                                                                                                                                   |
| UniProt ID         | <a href="#">S100B_HUMAN</a>                                                                                                                                             |
| Immunogen Region   | Met1-Glu92                                                                                                                                                              |
| Immunogen Sequence | QRFAQAQQQLPLESLGWDVA ELQLNHTGPQQDPRLYWQGG PALGRSFLHGPGLDKGQLRI<br>HRDGIYMVHIQVTLAICSST TASRHHPTTLAVGICSPASR SISLLRLSFHQGCTIASQRL<br>TPLARGDTLCTNLTGTLPS RNTDETFFGVQWVRP |

### ADDITIONAL INFORMATION

|      |                                                                                                                         |
|------|-------------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|-------------------------------------------------------------------------------------------------------------------------|



Recombinant Rat IL-4 Protein was determined by SDS-PAGE with Coomassie Blue, showing a band at 40-45kDa.