

## Human TACTILE/CD96 protein (Recombinant) (C-His) (STJP002919)

### GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

### PRODUCT PROPERTIES


Formulation	Lyophilised from a 0.2 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">10225</a>
Gene Symbol	<a href="#">CD96</a>
UniProt ID	<a href="#">TACT_HUMAN</a>
Immunogen Region	Val22-Met503
Immunogen Sequence	EEPLVVKVEEGDNAVLQCLK GTS DGPTQQLTWSRESPLKP FLKLSLGLPGLGIHMRPLAI WLFIFNVSQQMGGFYLCQPG PPSEKAWQPGWTVNVEGSGE LFRWNVSDLGGLGCLKNRS SEGPSSPSGKLMSPKLYVWA KDRPEIWEGEPCLPPRDSL NQSLSQDLTMAPGSTLWLSC GVPPDSVSRGPLSWTHVHPK GPKSLLSLELKDDRPARDMW VMETGLLLPRATAQDAGKY

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



Recombinant Human CD96 Protein was resolved with SDS PAGE under reducing (R) and non-reducing (NR) conditions, showing single bands at 90-140 kDa and 180-280 kDa, respectively.