

Human Complement C5A protein (Recombinant) (C-His) (STJP002921)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	E.coli

PRODUCT PROPERTIES


Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	727
Gene Symbol	C5
UniProt ID	CO5_HUMAN
Immunogen Region	Thr678-Arg751
Immunogen Sequence	ALPEDGGAAFPPGHFKDPKR LYCKNGGFFLRIHPDGRVDG VREKSDPHVKLQLQAEERGVS VSIKGVCANRYLAMKEDGRL LASKCVTEECFFFERLESNN YNTYRSRKYSSWYVALKRTG QYKLGSKTGPGQKAILFLPM SAKS

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



Recombinant Human Complement C5A Protein was determined by SDS-PAGE with Coomassie Blue, showing a band at 10-15kDa.