

Human EGLN1/PHD2 protein (Recombinant) (C-Flag) (STJP002923)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	E.coli

PRODUCT PROPERTIES


Formulation	Supplied as a 0.22 µg/ml filtered solution in 20mM Hepes, 250mM NaCl, pH 7.5.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	54583
Gene Symbol	EGLN1
UniProt ID	EGLN1_HUMAN
Immunogen Region	Gly177-Phe426
Immunogen Sequence	MLTELEKALNSIIDVYHKYS LIKGNFHAVYRDDLKLLKET ECPQYIRKKGADVWFKELDI NTDGAVNFQEFLILVIKMGV AAHKKSHEESHKE

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



Recombinant Human EGLN1/PHD2 Protein was determined by SDS-PAGE with Coomassie Blue, showing a band at 30-35kDa.