

Human MMP-7 protein (Recombinant) (C-His) (STJP002931)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES


Formulation	Supplied as a 0.22 µg/ml filtered solution in 10mM HEPES 150mM NaCl 10%Glycerol PH7.5.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	4316
Gene Symbol	MMP7
UniProt ID	MMP7_HUMAN
Immunogen Region	Tyr97-Lys267
Immunogen Sequence	PVYELQVQKSVTVQEGLCVL VPCSFSPWRSWYSSPLYV YWFRDGEIPYYAEVVATNNP DRRVKPETQGRFRLGDEVK KNCSLSIGDARMEDTGSYFF RVERGRDVKYSYQQNKLNLE VTALIEKPDHIFLEPLESGR PTRLSCLPGSCEAGPPLTF SWTGNALSPLDPETTRSEL TLTPRPEDHGTNLTCQMKRQ GAQVTTERTVQLNVSYAPQT ITIFRNGIALEILQNTSYL

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



Recombinant Human MMP7/MPSL1/PUMP1
 Protein was determined by SDS-PAGE with
 Coomassie Blue, showing a band at 35-40 kDa.