

Mouse CCL11/Eotaxin protein (Recombinant) (No-Tag) (STJP002939)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	E.coli

PRODUCT PROPERTIES


Formulation	Lyophilised from 0.22 Mu m filtered solution in PBS, 250mM NaCl, pH 7.4. Normally 8% Trehalose is added as protectant before lyophilisation.
Storage Instruction	Store at-20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	20292
Gene Symbol	Ccl11
UniProt ID	CCL11_MOUSE
Immunogen Region	His24-Pro97
Immunogen Sequence	TSHLVKCAEKEKTFCVNGGE CFMVKDLSNPSRYLCKCQPG FTGARCTENVPMKVQNQEKA EELYQKRVLTI

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



Recombinant Mouse CCL11/Eotaxin Protein was determined by SDS-PAGE with Coomassie Blue, showing a band at 10-15kDa.