

Mouse PPIase A/PPIA protein (Recombinant) (No-Tag) (STJP002942)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	E.coli

PRODUCT PROPERTIES


Formulation	Supplied as a 0.22 µg/ml filtered solution in PBS, 10% Glycerol, pH 7.5.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	268373
Gene Symbol	Ppia
UniProt ID	PPIA_MOUSE
Immunogen Region	Met1-Leu164
Immunogen Sequence	MDVFMKGLSKAKEGVVAAA E KTKQGVAAEAGKTKEGVLYV GSKTKEGVVHGVATVAEKT EQVTNVGGAVVTGVTAVAQK TVEGAGSIAAATGFVKKDQL GKNEEGAPQEGILEDMPVDP DNEAYEMPSEEGYQDYEEA

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



Recombinant Mouse PPIase A/PPIA Protein was determined by SDS-PAGE with Coomassie Blue, showing bands at 15-25 kDa.