

Human Alpha-synuclein/SNCA protein (Recombinant) (C-His) (STJP002951)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	Pichia

PRODUCT PROPERTIES


Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	6622
Gene Symbol	SNCA
UniProt ID	SYUA_HUMAN
Immunogen Region	Met1-Ala140
Immunogen Sequence	GPYGANVEDSICCQDYIRHP LPSRLVKEFFWTSKCRKPG VVLITVKNRDICADPRQVWV KKLLHKLS

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



Recombinant Human Alpha-synuclein/SNCA Protein was determined by SDS-PAGE with Coomassie Blue, showing a band at 19 kDa.