

Human S100-A10 protein (Recombinant) (C-His) (STJP002958)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES


Formulation	Supplied as a 0.22 µg/ml filtered solution in 20mM Tris 10% Glycerol, pH8.5
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	6281
Gene Symbol	S100A10
UniProt ID	S10AA_HUMAN
Immunogen Region	Met1-Lys97
Immunogen Sequence	GRETPDFGEVFDLDQQVWIF RNQALVTVP RSHRVTPVSVT ILPCKYPESLEQDKGIAIYL GIQNPDKCLFCKEVNGHPTL LLKEEKILDLYHHPEPMKPF LFYHTRTGGTSTFESVAFPG HYIASSKTGNPIFLT SKKGE YYNINFNLDIKS

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



Recombinant Human S100A10 Protein was determined by SDS-PAGE with Coomassie Blue, showing a band at 12.6 kDa.