

Human Integrin alpha V beta3 (ITGB3 & ITGAV) protein (Recombinant) (C-His) (STJP002962)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 M µm filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3685 3690
Gene Symbol	ITGAV ITGB3
UniProt ID	ITAV_HUMAN ITB3_HUMAN
Immunogen Region	Phe31-Val992
Immunogen Sequence	EAQTETCTVAPRERQNCGFP GVTPSQCANKGCCFDDTVRG VPWCFYPNTIDVPPEEECEF

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---