

Human IL-3R Beta/CSF2RB/CD131 protein (Recombinant) (No tag) (STJP002966)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	E.coli

PRODUCT PROPERTIES


Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1439
Gene Symbol	CSF2RB
UniProt ID	IL3RB_HUMAN
Immunogen Region	Trp17-Trp443
Immunogen Sequence	AELRCMCIKTTSGIHPKNIQ SLEVIGKGTHCNQVEVIATL KDGRKICLDPDAPRIKKIVQ KKLAGEDESAD

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



Recombinant Human IL-3R Beta/CSF2RB/CD131 Protein was determined by SDS-PAGE with Coomassie Blue, showing a band at 60-70 kDa.