

Human BCEI/PS2/TFF1 protein (Recombinant) (C-6His) (STJP002971)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	Pichia

PRODUCT PROPERTIES


Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	7031
Gene Symbol	TFF1
UniProt ID	TFF1_HUMAN
Immunogen Region	Glu25-Phe84
Immunogen Sequence	VVASELRCQCLKLPRVDFK NIQSLSVTPPGPHCAQTEVI ATKGGQKVCCLDPEAPLVQK IIQKILNKGKAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



Recombinant Human BCEI/PS2/TFF1 Protein was determined by SDS-PAGE with Coomassie Blue, showing a band at 17 kDa.