

Human HRG/HPRG protein (Recombinant) (C-His) (STJP002975)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	Pichia

PRODUCT PROPERTIES


Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3273
Gene Symbol	HRG
UniProt ID	HRG_HUMAN
Immunogen Region	Val19-Lys525
Immunogen Sequence	APMGSDPPTSCCFYSYTSRQL HRSFVMDYYETSSLCSKPAV VFLTKRGRQICANPSEPWVT EYMSDLELN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	-------------------------------------------------------------------------------------------------------------------------



Recombinant Human HRG/HPRG Protein was determined by SDS-PAGE with Coomassie Blue, showing a band at 75 kDa.