

Human CXCL7/PPBP protein (Recombinant) (C-hFc&6*His) (STJP002976)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES


Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5473
Gene Symbol	PPBP
UniProt ID	CXCL7_HUMAN
Immunogen Region	Ala59-Asp128
Immunogen Sequence	MDMFPLTWVFLALYFSRHQV RGQDPDPPCGGRLNSKDAGYI TSPGSPQDYP SHQNC EWIVY APEPNQKIVLNFNPHFEIEK HDCKYDFIEIRDGDSESADL LGKHCGNIAPPTIISSGSML YIKFTSDYARQGAGFSLRYE IFKTGSEDCSKNFTSPNGTI ESPGFPEKYPHNLDCTFTIL AKPKMEILQLIFDLEHDP LQVGEGDCKYDWLDIWDGIP HVGPLIGKYCGTKTPSELR

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



Recombinant Human CXCL7/PPBP Protein was determined by SDS-PAGE with Coomassie Blue, showing a band at 10-15 kDa.