

Human Ubiquitin protein (Recombinant) (C-His) (STJP002977)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES


Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	7314
Gene Symbol	UBB
UniProt ID	UBB_HUMAN
Immunogen Region	Met1-Gly76
Immunogen Sequence	ERVGCDLQPVDPTRGEVTF T SQVSEGCVAQAANAVREVH VLFLDFPGMLSHLELTQAS KQNGTETQEVFLVLSNKNV FVKFQAPEIPLHLAYDSSLV IFQGQPRVNITVLPSTSRK QILDWAATKGAITSIAALDD PQSIVLQLGQDPKAPFLCLP EAHKDMGATLEWQPRAQTPV QSCRLEGVSGHKEAYILRIL PGSEAGPRTVTVMMELSCTS GDAILILHGPPYVSWFIDI

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



Recombinant Human Ubiquitin protein was determined by SDS-PAGE under reducing conditions with Coomassie Blue.